



PRODUCT CATALOG

PAPTIZER SMART SANITIZER

SKU: T5L



Revolutionize your daily routine with the **LiViliti Paptizer Smart Sanitizer**, the ultimate sanitizing solution with cutting-edge technology. Designed to combine efficiency, safety, and convenience, the **Paptizer** utilizes UV-C light, heat, and circulated air to offer three distinct modes guaranteeing thorough sanitizing in as little as 3 minutes. Say goodbye to 99% of bacteria with just the press of a button, all without the use of ozone or mercury, making it the most advanced sanitizer on the market.

HOW IT WORKS

The Paptizer emits UV-C energy, absorbed by pathogens' nucleic acids, disrupting their reproduction and infectivity. This device utilizes peak germicidal wavelengths, rendering pathogens inactive for thorough disinfection

SPECIFICATIONS

PRODUCT DIMENSIONS: 11.89" L x 10.59" W x 6.06" H

TRAY DIMENSIONS: 9.45" L x 7.68" W x 4.72" H

PRODUCT WEIGHT: 4.29 LBS.

SANITIZING METHOD: UV-C LED

UV-C WAVELENGTH: 260-280 NM

UV-C LED QUANTITY: 40 PIECES

VOLTAGE: 110 VAC (60 HZ)

OPERATING TEMPERATURE: 32-104°F

HUMIDITY: 5-95%

WARRANTY: 2 YEARS



SCAN FOR VIDEOS



PAPTIZER SMART SANITIZER

KEY FEATURES

- **Ensures Cleanliness:** Personal care items can harbor bacteria and other harmful pathogens if not properly cleaned. The Paptizer's powerful disinfection capabilities eliminate over 99% of these microorganisms in just three minutes.
- **Promotes Health:** Clean equipment is essential for maintaining health and preventing infections. By effectively killing bacteria, the Paptizer helps reduce the risk of respiratory issues and other health problems associated with unclean equipment.
- **Easy to Use:** Users often face challenges when it comes to cleaning their equipment regularly. The Paptizer simplifies the sanitization process with its user-friendly operation, featuring three modes for sanitizing and simple one-touch buttons. This makes it easy for users to incorporate regular cleaning into their routine without hassle.
- **Comprehensive Coverage:** The Paptizer's multidirectional sanitization and 360-degree encircling certifies that every part of the equipment receives thorough coverage. This comprehensive coverage helps maintain cleanliness and hygiene, minimizing the risk of contamination.
- **Safe and Reliable:** Safety is a top priority and the Paptizer delivers with its ozone and mercury-free design. UV-C light, known for its safety and effectiveness in disinfection, provides a reliable solution for sanitizing equipment without compromising safety.
- **Long-lasting Performance:** With a 10,000+ hour LED lifespan and a 2-year warranty, the Paptizer is built to last, offering consumers reliable performance and peace of mind. Its durable design provides consistent sanitization for years to come, reducing the need for frequent replacements.
- **Portable Design:** The Paptizer's portable design allows users to take it wherever they go, ensuring access to clean and sanitized equipment even while traveling. Its compact size makes it convenient to use at home or on the go, providing users with flexibility and convenience.
- **Laboratory Tested:** Tested against common pathogens such as E. Coli, Staphylococcus aureus, Pseudomonas aeruginosa, and Candida albicans, the Paptizer delivers proven effectiveness in disinfection.